



## PROBES FOR LIQUID FEEDING



SKIOLD MAKES THE DIFFERENCE!

- This version is only used in long trough, where there is space for all the pigs at the same time.
- Pigs from 15 to 110 Kg is well served with this kind of probe.

Item No.	Name	Weight
140 309 050 790	Probe complete without mounting screws.	3,5 kg

4 sets of screws, bolt and washer is needed, but must be ordered separate depending on the type and thickness of the penning.

1. Mount the metal frame on the penning, with M8 Screws, of a good length, fit to thickness of the penning.

2. Mount the probe to the frame with the 2 clamps delivered together with the probe.


**Mounting:**  
Bolt+washer  
- penning plate -  
nut in front



The metal wire, is connected to the metal frame by just put it through the hole and then the bolt.

The Distance from the probe end to the bottom metal - default 25 mm.

Metal bottom part must be at the very bottom of the trough.



The wire from the frame - short to end of pipe. Connected with silicone connection. "Negative"

The wire from the probe end - 10 cm from pipe end. Connected with silicone connection. "Positive"

